

MIDMARK® 800 DIGITAL VITAL SIGNS DEVICE

	4-000-0801	4-000-0802	4-000-0803	4-000-0804
NIBP Mode: Spot	•	•	•	•
NIBP Mode: Average	•	•	•	•
NIBP Mode: SPRINT BP Protocol	•	•	•	•
NIBP Mode: Custom BP Protocol	•	•	•	•
Linear + Step Deflation Methods	•	•	•	•
Thermometry: Exergen® TemporalScanner™	•	•	•	•
Pain Score	•	•	•	•
Height	•	•	•	•
Pulse Rate	•	•	•	•
SpO ₂ Sensor: Masimo®			•	•
SpO ₂ Sensor: Nellcor®		•		
Automated Respiration Rate				•

Exergen, TemporalScanner, Masimo and Nellcor are trademarks of their respective owners.