



MEDICAL

Synthesis® Mobile Cabinets + Carts



Better care starts with a better design



Many healthcare organizations are realizing that better care starts with a better clinical environment design—one that keeps patients at the center of the care experience. We can help you take care to the patient, wherever that may be.



Make the most of the space you have

In a perfect world, we'd have all the space we need and everything within reach. **Synthesis® Mobile Treatment Cabinets by Midmark** bring the supplies you need to the point of care with specialized designs that help eliminate dead space and improve organization. And when you don't need them, the mobile cabinets can be easily stored out of the way.



STANDARD FEATURES



Concealed Hinges
Concealed, self-closing, nickel-plated metal hinges are standard with a clip-on mount, 110° opening and three adjustment points.



Seamless Drawers
Seamless, one-piece polystyrene drawers with rounded corners can help contain spills and be easily cleaned.



Casters
Choose from 60 mm casters and 101.6 mm quiet casters.



Seamless Panels
Seamless thermofoil panels are designed to resist spills, chemical degradation and delamination.




Drawer Slides
100-pound capacity, full extension slides can allow you to use all available drawer space.

Mobile Treatment Cabinet Options



Pulls
Choose from a variety of pulls including integrated, recessed and EPA-registered, antimicrobial options.
EPA registration number 85353-2



Security
Keyed and digital keyless lock options are available for additional security.



Top Styles Available
Kydex® Square Top
Kydex Compact Top
Kydex Contour Top



*Kydex Contour Top shown
(Available in .5" and 1.5" overhang)*

Soft Edge Bumper Top



Soft Edge Handle Top



Easily move from room to room

Care happens in many places within a healthcare facility. Caregivers are oftentimes required to move from room to room and patient to patient rapidly. To be effective, a mobile solution must allow for sturdy room-to-room maneuverability, even in difficult or small spaces. **Midmark® Mobile Procedure Carts** support procedure work and prep from the standing position with the flexibility in options to design the cart to work the way you need it.

STANDARD FEATURES



No-Tip Base

The no-tip base design reduces the likelihood of the cart tipping or falling.



Slide-Out Worksurface

The left or right slide-out worksurface provides an additional working area or space to place needed items.

Mobile Procedure Cart Options



01 Transit Handle
Provides a gripping surface for maneuvering the procedure cart



02 Soft Edge Contour Top with Raised Corners
Provides a barrier for the top surface and can be used to maneuver the cart



03 Sharps Bracket
Secure a sharps container to the procedure cart



04 Disinfectant Wipe Bracket
Secure a disinfectant wipe dispenser to the procedure cart (available in small and large sizes)



05 Waste Basket Bracket
Secure a waste basket to the procedure cart



06 6 x 6 Storage Bins
Add additional, removable storage to the procedure cart



07 Single Glove Holder
Secure a single box of sterile gloves to the procedure cart (top or side-mounted)



08 3-Box Glove Holder
Secure up to three boxes of sterile gloves to the procedure cart (top or side-mounted)



Kydex® Contour Top



Soft Edge Contour Top with Raised Corners



Keep essential tasks **moving** with care.

Cleaning, linen changes and waste removal are part of every busy exam and procedure room.

Midmark Synthesis® Mobile Linen and Waste Carts

help support efficient room-to-room workflows with versatile placement options—beneath or beside Synthesis Base Cabinets, or next to patient seating—blending seamlessly into the clinical environment.

When not in use, the carts store neatly out of the way and remain easily accessible, discreetly concealed with a clean, modern aesthetic that supports staff access and patient flow.

Mobile Linen Cart **Features**

Designed to streamline cleaning and linen removal for patient room beds and same-day procedure prep, the **Midmark Synthesis® Mobile Linen Carts** feature a hands-free, foot-activated lift for effortless linen removal and front lock for secure, sanitary transport in patient care areas. Locking casters ensure secure placement while two height options—31" and 33"—provide flexibility to fit a variety of clinical spaces. Front-opening design allows easy loading and unloading of linens for faster, more efficient room turnover.



Mobile Waste Cart **Features**

A cleaner-looking alternative to open trash solutions, the **Midmark Synthesis® Mobile Waste Carts** help keep waste concealed while supporting room-to-room use. Locking casters provide secure placement, while two height options—31" and 33"—offer flexibility across a variety of clinical spaces. Front-access waste design enables quick removal and replacement to support clean, efficient workflows.



Choose a style that complements your unique brand.

The Synthesis® Cabinetry Collection by Midmark includes many styling options including colors, finishes and other elements designed to allow you to put your unique signature on a space, whether one room or a network of facilities, for a consistent brand perception.



COLOR OPTIONS

	Mist 980*
	Venus Silver 804
	Flat Iron 807
	Storm 802
	Riverstone 967
	Pearl Essence 805
	Ember 968
	Tundra 966
	Frost 791, 931*
	Fog 969*
	Pebble Grey 216, 793*
	Sandbar 963
	Prairie 622
	Dune 794
	Flax 975
	Oat 611
	Russet 964
	Timber 801
	Sable 965

*PVC Free available.

HANDLE STYLES



**EPA registration number 85353-2



Designing better care.®



USGBC
MEMBER

CARB 93129.2 Phase 2 Compliant and TSCA Title VI Compliant

No product can guarantee points. Points are cumulative for all effected materials for a LEED® project. For further information on LEED, please visit the USGBC website: usgbc.org. LEED is a registered trademark of the US Green Building Council.

Kydex is a registered trademark of SEKISUI POLYMER INNOVATIONS, LLC

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: midmark.com/ISO

For more information, contact your Midmark dealer or call: 1.800.MIDMARK Outside the USA call: 1.937.526.3662 or visit our website: midmark.com

© 2026 Midmark Corporation, Versailles, Ohio USA
Products subject to improvement changes without notice.
007-10288-00 Rev. D1 (6/26)

