

Mobile Lift Table

Lift.
Prep.
Transport.



The Midmark Mobile Lift Table can be placed where you need it most, offering the utmost flexibility within your facility. Ergonomically designed to support safety and care team convenience, it's an ideal table for a variety of clinical care situations.

- Prevent fluids from running off the table with a raised marine edge.
- Protect staff hips and ankles with rounded corners on tabletop and base.
- Lift tabletop to proper ergonomic position with quiet, hands-free foot switch.
- Enhance efficiency with a battery located under the tabletop for improved connectivity and accessibility for cleaning.
- Enjoy greater stability with enhanced tube-frame table construction.
- Protect walls with non-marring bumpers.
- Move table easily with corrosion-resistant nylon locking casters
- Restrain restless patients with tie-down cleats.
- Enhance workflow with integrated scale.



MIDMARK MOBILE
LIFT TABLE SPECIFICATIONS



- Length:** 44.5"
- Width:** Tabletop 20" and Base frame 30.75"
- Height:** Without Scale - Adjustable from a low of 14.25" to a maximum of 36.25"
- With Scale - Adjustable from a low of 16.5" to a maximum of 38.5"
- Capacity:** Lifts up to 200 lb
- Battery:** Enhanced battery capacity holds charge for up to two weeks with normal use.
- Optional:** Integrated Scale

ORDERING GUIDE

Part Number	Description:
903-3120-04	Mobile Lift Table without Scale
903-3120-05	Mobile Lift Table with Scale

