



MONITOR MOUNT

IMPROVED EFFICIENCY. BETTER ERGONOMICS. MORE WORKSPACE.

The Midmark sliding track monitor mount is designed to assist dental clinicians with positioning a monitor on a Synthesis® and Artizan® Expressions rear cabinet, providing an efficient and ergonomic method of viewing media such as office schedules, patient charts and digital x-rays. The mount allows the monitor to be suspended above the countertop, providing more space to work and extra room for essential supplies and ancillary devices.

FEATURES



SIDE-TO-SIDE AND WRAP-AROUND

Sliding track allows mount to slide from doctor's side to assistant's side and also wrap around the cabinet to be viewed outside the operatory.



PIVOT

Left-right pivoting of monitor assists dental clinicians in viewing media from a variety of angles.



TILT

Tilt feature helps reduce glare when viewing from a standing or seated position.



HEIGHT ADJUSTMENT

Height adjustment mechanism allows the monitor to be adjusted to avoid support trays and taller devices.

ADDITIONAL FEATURES

- Easy positioning of monitor
- Accommodates both 100/75 VESA mounting patterns
- Finished back cover
- Cable management system
- Cable options available

SPECIFICATIONS

- Maximum monitor weight: 15 lb
- Accommodates up to 19" monitor
- Vertical adjustment: 5"
- Swivel range: 180°
- Tilt range: 15° forward, -15° backward
- Track length: 22"

**2-YEAR
LIMITED
WARRANTY**