




MEDICAL

Synthesis® Mobile Cabinets + Carts



Better care starts with a better design



Many healthcare organizations are realizing that better care starts with a better clinical environment design—one that keeps patients at the center of the care experience. We can help you take care to the patient, wherever that may be.



Make the most of the space you have

In a perfect world, we'd have all the space we need and everything within reach. **Synthesis® Mobile Treatment Cabinets by Midmark** bring the supplies you need to the point of care with specialized designs that help eliminate dead space and improve organization. And when you don't need them, the mobile cabinets can be easily stored out of the way.



STANDARD FEATURES



Concealed Hinges
Concealed, self-closing, nickel-plated metal hinges are standard with a clip-on mount, 110° opening and three adjustment points.



Seamless Drawers
Seamless, one-piece polystyrene drawers with rounded corners can help contain spills and be easily cleaned.



Casters
Choose from 60 mm casters and 101.6 mm quiet casters.



Seamless Panels
Seamless thermofoil panels are designed to resist spills, chemical degradation and delamination.



Drawer Slides
100-pound capacity, full extension slides can allow you to use all available drawer space.

Mobile Treatment Cabinet Options



Pulls
Choose from a variety of pulls including integrated, recessed and EPA-registered, antimicrobial options.
EPA registration number 85353-2



Kydex® Contour Top



Soft Edge Bumper Top



Security
Keyed and digital keyless lock options are available for additional security.



Soft Edge Handle Top

Easily move from room to room

Care happens in many places within a healthcare facility. Caregivers are oftentimes required to move from room to room and patient to patient rapidly. To be effective, a mobile solution must allow for sturdy room-to-room maneuverability, even in difficult or small spaces. **Midmark® Mobile Procedure Carts** support procedure work and prep from the standing position with the flexibility in options to design the cart to work the way you need it.

STANDARD FEATURES



No-Tip Base

The no-tip base design reduces the likelihood of the cart tipping or falling.



Slide-Out Worksurface

The left or right slide-out worksurface provides an additional working area or space to place needed items.

Mobile Procedure Cart Options



01 Transit Handle
Provides a gripping surface for maneuvering the procedure cart



02 Soft Edge Contour Top with Raised Corners
Provides a barrier for the top surface and can be used to maneuver the cart



03 Sharps Bracket
Secure a sharps container to the procedure cart



04 Disinfectant Wipe Bracket
Secure a disinfectant wipe dispenser to the procedure cart (available in small and large sizes)



05 Waste Basket Bracket
Secure a waste basket to the procedure cart




06 6 x 6 Storage Bins
Add additional, removable storage to the procedure cart



07 Single Glove Holder
Secure a single box of sterile gloves to the procedure cart (top or side-mounted)



08 3-Box Glove Holder
Secure up to three boxes of sterile gloves to the procedure cart (top or side-mounted)



Kydex® Contour Top



Soft Edge Contour Top with Raised Corners





Designing better care.®



CARB 93129.2 Phase 2 Compliant and TSCA Title VI Compliant

No product can guarantee points. Points are cumulative for all effected materials for a LEED® project. For further information on LEED, please visit the USGBC website: usgbc.org. LEED is a registered trademark of the US Green Building Council.

Kydex is a registered trademark of SEKISUI POLYMER INNOVATIONS, LLC

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: midmark.com/ISO

For more information, contact your Midmark dealer or call: 1.800.MIDMARK Outside the USA call: 1.937.526.3662 or visit our website: midmark.com

© 2020 Midmark Corporation, Miamisburg, Ohio USA
Products subject to improvement changes without notice.
Litho in USA.
007-10288-00 Rev. C1 (7/22)

