




MEDICAL

# Digital Diagnostics Connected Solutions



# You don't have to settle for complexity.




At Midmark, we believe technology should work the way you do. That's why we designed our digital diagnostic solutions to simplify the complex—bringing devices, software and support together in one connected platform.



# Designed with connectivity in mind

Your Midmark vital signs device is designed for seamless connection to your EMR platform.



Midmark Connect is the bridge between your diagnostic devices and the EMR—developed to help eliminate manual entry, streamline clinical workflows and support more accurate, timely patient data.

With Midmark Connect, diagnostic data flows the way it should—securely, accurately and without workarounds—so your team can focus more on patients, not technology.



SCAN TO SEE  
THE WORKFLOW

# The Midmark Vital Signs Workflow




### LAUNCH THE TEST

01 Open the EMR from the patient record and select launch.



### ACQUIRE THE TEST

02 Acquire the diagnostic test without leaving the EMR.



### REVIEW + SAVE

03 Verify the test results and save.



### SYNC DATA TO THE EMR

04 Test data transfers directly into the EMR.

*Workflow may vary slightly among EMRs.*

## SEAMLESS COMPATIBILITY. CONFIDENCE FROM THE START.

Disconnected solutions force you to jump between systems, manage middleware or re-enter results by hand. That slows care, frustrates staff and increases the risk of errors. Midmark helps eliminate those barriers.

- Never leave the patient chart
- No middleware required
- Streamlined user workflow across vital signs and ECG
- Seamless transcription of results to the EMR
- Already configured to integrate with today's leading EMRs



# We've thought through the risks—so you can focus on care.

Bringing new diagnostic devices into your practice can feel risky. Will they connect with your EMR? Will they meet your IT and compliance standards? Will they truly work the way you do? With Midmark, you don't have to wonder. Our solutions are built to integrate with the EMRs you already know and trust. If your practice isn't presently using an EMR, Midmark IQmanager® can bridge the gap to keep your workflow connected.

## Product Cybersecurity by Design

Cybersecurity is integrated throughout the Midmark product lifecycle to help protect device functionality and patient data. Through risk-based design, security testing, vulnerability management and ongoing support, we help healthcare organizations securely deploy and maintain connected Midmark solutions.

## Proof of Concept

And before you make any commitment, we'll help you test and validate our solutions in your own environment. Through our Proof of Concept program, you can confirm integration, performance and fit with your clinical and IT teams—reducing risk and giving you confidence before you move forward. Connect with your Midmark Sales Representative to see if you qualify.



SCAN TO SEE  
COMPATIBLE EMRS

NIST Cybersecurity Framework is referenced for alignment only.

# The Right Device for Your Unique Workflow

Midmark offers a portfolio of digital diagnostic devices—built for accuracy, efficiency and seamless EMR connectivity.



**Midmark® 800 Digital Vital Signs Device**

- Automates blood pressure, temperature and pulse rate with optional SpO<sub>2</sub> and respiration rate. Connect the Midmark 800 to the Midmark 626 Barrier-Free® Examination Chair with an integrated digital scale to capture weight.
- Linear + Step Deflation NIBP methods.
- Built in SPRINT protocol with patient positioning prompts as well as spot, custom and averaging BP modes
- Large, high-contrast screen for easy viewing.
- Operate a connected exam chair and digital scale from a single screen — the device or a computer.
- Brings automation to the point of care—designed for consistent, repeatable vital signs measurements.
- Seamlessly sends data to the EMR, reducing the likelihood of manual transcription errors.



**Midmark® Digital ECG (with optional DICOM compatibility)**

- Simple acquisition with automatic EMR integration
- DICOM compatibility supports secure data exchange between the care team's EMR and the cardiologist's preferred CVIS
- May help reduce manual transcription and speeds diagnostic decision-making
- Designed for accuracy and efficiency in daily use

**Whichever solution you choose, Midmark devices work the way you do—streamlining workflows while supporting accurate, reliable patient data capture.**

# From device to delivery, Midmark simplifies purchasing.

Every environment is unique—which is why we offer flexible mounting and purchasing options tailored to your workflow.

**MOUNTING SOLUTIONS FOR VITAL SIGNS**

**MIDMARK® 800 DIGITAL VITAL SIGNS DEVICE**



Mobile Cart  
PN: 9A747001



Wall Board Mount  
PN: 9A748001



Mobile Workstation Mount  
PN: 002-11619-00

**PURCHASING PROCESS MADE SIMPLE**

Hardware (devices, carts, mounts) is purchased through authorized Midmark distributors. Software (licensing, integration, implementation) is purchased directly from Midmark and requires a quote from your Midmark Sales Representative. *(Mounting and software sold separately.)*

With our coordinated process, hardware and software orders are aligned to support a smooth installation and fast go-live.

**WHAT YOU’LL NEED:**

- Completed IT questionnaire
- Product and software configuration
- Site details for delivery and install planning



Because  
true  
partnership  
means  
standing  
with you  
long after  
the sale.

“Our commitment doesn’t end at installation. With Midmark, you get the confidence of knowing we’ll be there—through setup, training and the long life of your devices.”

#### MAINTENANCE + SUPPORT

Your first year of software Maintenance and Support is included with your quote. Annual renewals keep you secure, compliant and up to date with the latest features and patches. Midmark-trained technicians and reliable service programs help ensure your investment continues to perform day in and day out.

#### POST-SALE SUPPORT

From ongoing software updates and warranty coverage to clinical education and Live Design services, we’re here to help your team succeed. Whether it’s troubleshooting connectivity, planning your next space or training staff on best practices, you’ll always have access to the expertise and resources you need.

---

## SERVICES + SUPPORT

#### CLINICAL EDUCATION

Trust in our expert educators to assist you in implementing your Midmark devices.



SCAN FOR  
TRAINING  
OPTIONS.

#### MIDMARK LIVE DESIGN

Our in-house design experts are ready to help you every step of the way.



SCAN TO  
VIEW VIDEO.

#### REPAIR + SERVICE

Explore available support and service plans so you can obtain the level of coverage you need to keep your medical equipment and devices operational. This includes periodic maintenance and periodic inspection plans, as well as extended warranty plans.



SCAN FOR  
SUPPORT.

#### CONTACT US

Have questions or ready to get started? Connect with a Midmark representative for software licensing information, detailed system requirements and to see what software you will need.



SCAN TO  
CONNECT.



Designing better care.®

CARB 93120.2 Phase 2 Compliant and TSCA Title VI Compliant

Midmark is an ISO 13485 and ISO 9001 Certified Company.  
Certain products are not included. See the complete list at:  
[midmark.com/ISO](http://midmark.com/ISO)

For more information, contact your authorized Midmark  
Distributor or call: 1.800.MIDMARK

Outside the USA call: 1.937.526.3662 or visit our website:  
[midmark.com](http://midmark.com)

© 2026 Midmark Corporation, Versailles, Ohio USA  
Products subject to improvement changes without notice.

007-10543-00 Rev. B1 (8/26)

