



FURNITURE IDIQ GOVT CONTRACT

2025 Guide Book

General Information

Ordering Notes.....	3
Freight and Shipping.....	3

Product Information

MTR CODE	Series
CM3	Midmark Mobile Treatment Cabinets and Ritter Waste Receptacle
CM3	P Series Procedure Carts
CS4	6281 Laptop Wall Mount Workstation
CS4	6282 PC Wall Mount Workstation
CS4	6283 Heavy Duty PC Wall Mount Workstation
GC1	Midmark Seating 280 and 282 Chairs
GC2	Ritter 282 Bariatric Side Chairs
HAW4	Midmark Workstations
MM1	Midmark Instrument Panels, Dressing Nooks, Glove/Towel/Cup Dispensers, and Exam Room Cabinetry Package
MSS1	Midmark Synthesis Casework (layout designed to request)
MSS2	Midmark Synthesis Casework (layout designed to request)
PC1	280/282 Patient Side Chairs
PC3	282 Bariatric Chair
PP1	Midmark Synthesis Casework
SC1	Midmark Synthesis Casework
SH1	Midmark Synthesis Casework
ST1	425/427 Provider & Lab Stools
ST1	Ritter 27X-002 Stools without back
ST2	425/427 Provider & Lab Stools
ST2	Midmark 270 Series Stools w/Lab Extension & Foot Ring

Warranty Information

Midmark® and Ritter® Medical Equipment Warranty	36
Midmark Workstation Warranty.....	37
Midmark Synthesis Cabinetry Warranty.....	38

Midmark Customer Experience

For order status, requests and product information, contact Medical Customer Experience: **1.844.856.1233**

Email: **allcustservice-medical@midmark.com**

To place an order, email: **Orders-Medical@midmark.com**

For Synthesis® Cabinetry order placement or information, email: **cabinetry-medical@midmark.com**

Midmark Technical Service

For service part identification, technical/operational questions, and warranty and product issues/complaints, contact Technical Service: **1.844.856.1230**

Sterilizer + Ultrasonic Products: **sterilizers.ts@midmark.com**

Diagnostic Products: **techsupport@midmark.com**

Procedure + Exam Chairs, Seating, Lighting, Workstations, Diagnostic Instrument Wall/Desk Kits, Oscopes,

Ophthalmoscopes + Sphygmomanometers: **examroomeq.ts@midmark.com**

Special Notes

1. ORDERING—Please provide the following information on all orders:

- Purchase order number
- Bill-to and ship-to addresses
- Product models and/or part number(s)
- Description including upholstery color and/or panel/countertop color
- Quantity needed
- Pricing
- Required shipping date
- Special shipping requirements (see details below)

DROP-SHIP REQUIREMENTS

- Please specify if the customer has a loading dock or has special shipping requirements.
- Receiving contact name, email and phone number is required.

2. SPECIAL DESIGN—Orders for equipment of special design, construction, upholstery or special finish are not subject to cancellation after your acceptance of our acknowledgment.

3. SPECIAL SHIPPING REQUIREMENTS—These services are available to enhance delivery to your customers. Please call for pricing.

APPOINTMENT REQUIRED—Carrier will call delivery contact and mutually agree upon delivery appointment. This is not a guaranteed time/date.

LIFTGATE—Carrier uses liftgate/ramp to get product to ground level when facility does not have a dock. Facility is responsible for taking product once at ground level.

INSIDE DELIVERY—Carrier unloads product and delivers inside the first set of doors only, provided entrance to facility can accommodate freight as packaged (including skids) and there are no steps.

MIDMARK DELIVERY SERVICE—Carrier sets product in place per the delivery contact. Carrier is responsible for assembling product and removal of all packaging, trash and skids from the delivery site (Synthesis® Cabinetry does not include installation)

4. FREIGHT DAMAGE—Merchandise arriving with external damage should be so described on the bill of lading. In case of concealed damage, carrier must be notified within five days after receipt and claims for damage instituted.

5. ORDER CHANGES—All changes must be submitted in writing. Any changes made within five days of the order shipping will result in the order being moved to the next available ship date.

6. ORDER CANCELLATION—Any order, if cancelled within 10 days before the scheduled shipping date, will be charged a 20% cancellation charge.

7. RETURN OF MERCHANDISE—A Return Material Authorization (RMA) must be obtained for the return of any item. If return is authorized, transportation charges must be prepaid. Risk of loss with respect to returned goods shall remain with purchaser until receipt by Midmark. Items returned for credit are subject to a charge of 25% or more to cover inspection, damage in transit, reconditioning and repacking. Items may be returned within 6 months (except service parts within 90 days) from the original ship date if properly packaged and in a good, unused condition. RMA # should be written on the outside of the packaging. Diagnostic products may be returned within 30 days if in a good, unused condition and in original sealed packaging.

Note: Certain items are non-returnable: Specials, Synthesis Cabinetry, Midmark® Workstations, Discontinued Items or Parts, Custom or Non-Standard Colors, Service Upholstery and Service Parts if value is less than \$50, and Diagnostic Consumables; e.g., Electrodes, all Disposable Items and Diagnostic Items valued less than \$50 per item.


8. SPECIAL ORDER REQUIREMENTS—Any order requiring any special deviation from standard products will carry the following disclaimer:

We hereby accept and acknowledge the above captioned order. However, because of the special nature of this order in that this equipment is being supplied upon special design specifications, special materials and/or field modifications, as per your request, the standard warranties of Midmark and Ritter on its products are not applicable to this order, and this acceptance and acknowledgment is conditioned upon your indemnification and promise to hold Midmark and Ritter, its officers, employees and agents, harmless in respect to any and all warranty claims, damages, both direct and consequential, expenses, costs and any liability arising from the use of or relating to this special equipment.

9. IMPROVE PRODUCT QUALITY—Midmark Corporation reserves the right to change design and construction without notice.

Upholstery Color Palette

XXX - XXX - XXX























Example:
630-010 *Chair*
002-10108-843 *Top*
630 HUMANFORM® *Procedure Chair - Oasis, Premium Upholstery Top*

Example: 276-001-814
 276 *Physician Seating, Stone Upholstery*











Standard Upholstery Colors

Color is an important part of establishing your distinct brand. It's more than visual appeal. Color can transform your space into something you won't just see, but you will feel, too. We offer many upholstery color selections to complement the look, feel and intended function of your space.

	Mist-813		Healing Waters-848
	Stone-814		Shaded Garden-853
	Restful Path-858		Citrus-849
	Iron Ore-851		Sandy Retreat-850
	Lunar Gray-845		Curative Copper-846
	Obsidian-857		Cranberry-859
	Oasis-843		Dark Linen-856
	Soothing Blue-855		Latte-860
	Dream-847		Robust Brown-852
	Mineral-844		Deep Earth-854

UltraFree Upholstery Colors

There is nothing quite like the rich feel of leather, but it's too difficult to clean and maintain to be practical in the clinical environment. With UltraFree upholstery, you can have the look and feel of leather with PVC-free material. And now with seamless application, it is easier to clean and disinfect.

	UltraFree Stone-865		UltraFree Spice-863
	UltraFree Restful Path-868		UltraFree Cranberry-869
	UltraFree Obsidian-867		UltraFree Dark Linen-866
	UltraFree Harbor-861		UltraFree Latte-870
	UltraFree Wheatgrass-862		UltraFree Branch-864



Midmark® 425

Physician Seating (Hand Operated)

- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description

Basic stool w/o seat cushion
Seat cushion

Model

425-001
002-0273-XXX^{1,2}

Accessories

002-0274-XXX^{1,2}
002-0275-XXX^{1,2}
002-1546-00

Back kit
Arm kit
Lab height w/ black foot ring (425 only)

¹ - See color palette on page 4

² - Additional charge for Special or UltraFree Upholstery/Accessory Color

MTR CODE

ST1
ST2



Midmark® 427

Physician Seating (Foot Operated)

- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description

Basic stool w/o seat cushion
Seat cushion

Model

427-001
002-0273-XXX^{1,2}

Accessories

002-0274-XXX^{1,2}
002-0275-XXX^{1,2}

Back kit
Arm kit

¹ - See color palette on page 4

² - Additional charge for Special or UltraFree Upholstery/Accessory color

MTR CODE

ST1
ST2



Model

279-001-XXX^{1,2}

Ritter® 279

Physician Seating (Hand Operated)

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 30 lb
- Shipping weight: 35 lb

Description

Basic stool w/hand release - pneumatically adjustable

MTR CODE

ST2



Model

278-001-XXX^{1,2}

Ritter® 278

Physician Seating (Foot Operated)

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 30 lb
- Shipping weight: 35 lb

Description

Basic stool w/foot release - pneumatically adjustable

MTR CODE

ST2



hand release

Model

277-001-XXX^{1,2}
277-201-XXX^{1,2}

Ritter® 277

Physician Seating with Backrest

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description

Basic stool w/back - pneumatically adjustable, hand release
Basic stool w/back - pneumatically adjustable, foot release



foot release

MTR CODE

ST2



hand release

Model

276-001-XXX^{1,2}
276-201-XXX^{1,2}

Ritter® 276

Physician Seating

- Chrome five leg base
- Weight capacity: 350 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Description

Basic stool - pneumatically adjustable, hand release
Basic stool - pneumatically adjustable, foot release



foot release

MTR CODE

ST2

¹ - See color palette on page 4

² - Additional charge for Special or UltraFree Upholstery color



Model
275-001-XXX^{1,2}

Ritter® 275

Physician Seating with Backrest

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description

Basic stool w/back - manually adjustable

MTR CODE

ST2



Model
274-001-XXX^{1,2}

Ritter® 274

Physician Seating

- Chrome five leg base
- Weight capacity: 350 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Description

Basic stool - manually adjustable

MTR CODE

ST2



Model
273-001-XXX^{1,2}

Ritter® 273

Physician Seating with Backrest

- Black composite five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description

Basic stool w/back - pneumatically adjustable

MTR CODE

ST2



Model
272-001-XXX^{1,2}

Ritter® 272

Physician Seating

- Black composite five leg base
- Weight capacity: 300 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Description

Basic stool - pneumatically adjustable

MTR CODE

ST2

¹ - See color palette on page 4

² - Additional charge for Special or UltraFree Upholstery color



Model
271-001-XXX^{1,2}

Ritter® 271

Physician Seating with Backrest

- Black composite five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description
Basic stool w/back - manually adjustable

MTR CODE
ST2



Model
270-001-XXX^{1,2}

Ritter® 270

Physician Seating

- Black composite five leg base
- Weight capacity: 300 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Description
Basic Stool - manually adjustable

MTR CODE
ST2

Model
Miscellaneous

- 9A366001-XXX^{1,2}
- 27X-002-XXX^{1,2}
- 27X-003-XXX^{1,2}
- 27X-004-XXX^{1,2}
- 27X-005-XXX^{1,2}
- 9A367001
- 9A368001

Description

Back rest kit (weight capacity, 225 lb)

Set of soft rubber casters

Set of glides (not available on foot release models)

Set of locking casters

Auto locking casters

Foot ring, field installed (not available on foot release models)

6" lab height ext. kit (not available on foot release models)

MTR CODE
ST2
ST1/ST2
ST2
ST2
ST2
ST2
ST2

¹ - See color palette on page 4

² - Additional charge for Special or UltraFree Upholstery/Accessory color



Ritter® 280

Side Chair

- 400 lb patient weight capacity
- Wall-saver, kickback rear leg design
- Leveling feet
- Chair weight: 25 lb
- Shipping weight: 35 lb

MTR CODE

GC1
PC1

Model

Description

280-001-XXX ¹	Basic chair (without arms)
280-001-XXX ¹	Basic chair (without arms) UltraFree upholstery
280-002-XXX ¹	Basic chair (with arms)
280-002-XXX ¹	Basic chair (with arms) UltraFree upholstery

Miscellaneous

280-001-999	Basic chair Special Color (without arms)
280-002-999	Basic chair Special Color (with arms)

¹ - See color palette on page 4



Ritter® 282

Bariatric Side Chair

- 650 lb patient weight capacity
- Wall-saver, kickback rear leg design
- Leveling feet
- Chair weight: 37 lb
- Shipping weight: 47 lb

MTR CODE

GC1
GC2
PC1
PC3

Model

Description

282-001-XXX ¹	Basic chair
282-001-XXX ¹	Basic chair UltraFree upholstery

Miscellaneous

282-001-999	Basic chair Special Color
-------------	---------------------------

¹ - See color palette on page 4



Instrument Panel

- Instruments not included
- Shipping weight: 17 lb

MTR CODE
MM1

Product Information

Step 1. Configuration

- I1** Instrument Panel - 30" x 20"
- I2** Instrument Panel - 34" x 20"

Step 2. Option Package

- X** No Lock

Step 3. Cabinetry Finish*

- XXX** Solid Color Cabinet Finish (Excludes PVC-Free)
- XXX** Woodgrain, Metallic + PVC-Free Cabinet Finish

Configuration Example: I2X.803

*See page 24 for Cabinetry Finish

MTR CODE

HAW4



Midmark® 6211

Secure Laptop Workstation

- 26" of height adjustment from 23" to 49"
- Keyed laptop security system
- 23" D x 26" W worksurface
- 11" D x 20.5" W secondary worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- 19.5" D x 20.5" W base dimension
- Weight capacity: 20 lb
- Unit weight: 58 lb
- Shipping weight: 103 lb

Model

6211-001-xxx¹

Description

Secure Laptop Workstation

Accessories

Plant Installed

Field Installed

BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/ Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-00 ²	Equipment Pole Mount for Digital Vital Signs Device
MPSF	9A648001	Mobile Printer Shelf
PMSL ³	9A655001 ³	Power Mount w/ Steel Cable Mgmnt-LH
PMSR ³	9A655001 ³	Power Mount w/ Steel Cable Mgmnt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Power Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Power Cord-RH
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB ⁴	9A656001 ⁴	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount

¹ See color code chart on page 24

² Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.

³ 9A655001 & 9A655002 - Extended Rear legs required p/n 020-10016-00 for models 6211, 6214, 6215-001, 6215-002, 6217

⁴ 9A656001 - Extended Rear legs required p/n 020-10016-00 for models 6211, 6214, 6215-001, 6215-002, 6215-003, 6217



Midmark® 6212

Flat Panel PC Workstation

- 26" of height adjustment from 23" to 49"
- 23" D x 26" W worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical adj. and up to 40° positive tilt
- Monitor arm capacity range 5 lb min - 20 lb max
- Monitor arm independent rotation of 180°
- 22.5" D x 20.5" W base dimension
- Weight capacity: 30 lb
- Unit weight: 55 lb
- Shipping weight: 100 lb

Model

6212-001-xxx¹

Description

Flat Panel PC Workstation

Accessories

Plant Installed

Field Installed

AIOM	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLW	9A647001	Drawer w/ Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-002	Equipment Pole Mount for Digital Vital Signs Device
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
PMSL ³	9A655001 ³	Power Mount w/ Steel Cable Mgmt-LH
PMSR ³	9A655001 ³	Power Mount w/ Steel Cable Mgmt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Power Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Power Cord-RH
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB ⁴	9A656001 ⁴	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount

¹ See color code chart on page 24

² Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.



Midmark® 6213

Flat Panel Secure PC Workstation

- 26" of height adjustment from 23" to 49"
- Keyed PC security system
- 23" D x 26" W worksurface
- 11" D x 20.5" W secondary worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical adjustment and up to 40° positive tilt
- Monitor arm capacity range 5 lb min - 20 lb max
- Monitor arm independent rotation of 180°
- 22.5" D x 20.5" W base dimension
- Weight capacity: 30 lb
- Unit weight: 64 lb
- Shipping weight: 109 lb

Model

6213-001-xxx¹

Description

Flat Panel Secure PC Workstation

Accessories

Plant Installed

AIOM

BCHD

BCHL

BCHR

BPCH

CPUB

CPUF

DBCM

DMMT

DWLW

DWNW

EPWS

IQEP²

LPDH

MPSF

PCML

PCMR

PMBL³PMBR³

RBWS

SBCB

SBCF

SHLH

SHRH

TBDH

USCB⁴

VESA

VIPL

VIPR

VPBL

VPBR

Field Installed

N/A

9A664001

9A659001

9A659001

9A658001

9A663001

9A663001

9A653001

9A667001

9A647001

9A646001

9A651001

002-10592-00²

9A662001

9A648001

9A669001

9A669001

9A655002³9A655002³

9A650001

9A652001

9A652001

9A654001

9A654001

9A661001

9A656001⁴

9A660001

9A649001

9A649002

9A695001

9A695002

Monitor over 10 lb

Bar Code Holder Mount

Bar Code Holder-LH

Bar Code Holder-RH

Blood Pressure Cuff Holder

CPU Holder-Back

CPU Holder-Front

Double Basket

Dual Monitor Mount

Drawer w/ Lock-Work Surf Mnt

Drawer-Work Surf Mnt

Equipment Pole

Equipment Pole Mount for Digital Vital Signs Device

Laptop Display Holder

Mobile Printer Shelf

MG Coiled Power Cord + Steel Hanger - LH

MG Coiled Power Cord + Steel Hanger - RH

Power Mount w/ Bin Cable Mgmnt-LH

Power Mount w/ Bin Cable Mgmnt-RH

Rear Basket

Single Basket-Back

Single Basket-Front

Scale Holder-LH

Scale Holder-RH

Tablet Display Holder

Adjustable Utility Shelf

VESA Mounting Plate

VOIP Shelf-LH Mount

VOIP Shelf-RH Mount

VOIP Steel Bracket-LH Mount

VOIP Steel Bracket-RH Mount

¹ See color code chart on page 24² Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.

MTR CODE

HAW4



Midmark® 6214

Procedure Workstation

- 26" of height adjustment from 23" to 49"
- 21" D x 26" W worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Includes: rear basket, BP cuff holder, integrated spirometer holsters and equipment pole
- 19.5" D x 20.5" W base dimension
- Weight capacity: 30 lb
- Unit weight: 56 lb
- Shipping weight: 101 lb

Model

6214-001-xxx¹

Description

Procedure Workstation

Accessories

Plant Installed

Field Installed

CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/ Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
IQEP	002-10592-00	Equipment Pole Mount for Digital Vital Signs Device
PMSL ²	9A655001 ²	Power Mount w/ Steel Cable Mgmnt-LH
PMSR ²	9A655001 ²	Power Mount w/ Steel Cable Mgmnt-RH
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB ³	9A656001 ³	Adjustable Utility Shelf

¹ See color code chart on page 24
² 9A655001 & 9A655002 - Extended Rear legs required p/n 020-10016-00 for models 6211, 6214, 6215-001, 6215-002, 6217
³ 9A656001 - Extended Rear legs required p/n 020-10016-00 for models 6211, 6214, 6215-001, 6215-002, 6215-003 , 6217

MTR CODE

HAW4



Midmark® 6215-001

Rectangular Workstation

- 26" of height adjustment from 23" to 49"
- 21" D x 26" W worksurface
- 19.5" D x 20.5" W base dimension
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Weight capacity: 20 lb
- Unit weight: 48 lb
- Shipping weight: 93 lb

Model

6215-001-xxx¹

Description

Compact Workstation

Accessories

Plant Installed

Field Installed

BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/ Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-00 ²	Equipment Pole Mount for Digital Vital Signs Device
MPSF	9A648001	Mobile Printer Shelf
PMSL ³	9A655001 ³	Power Mount w/ Steel Cable Mgmt-LH
PMSR ³	9A655001 ³	Power Mount w/ Steel Cable Mgmt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Power Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Power Cord-RH
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB ⁴	9A656001 ⁴	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount

¹ See color code chart on page 24

² Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.

³ 9A655001 & 9A655002 - Extended Rear legs required p/n 020-10016-00 for models 6211, 6214, 6215-001, 6215-002, 6217

⁴ 9A656001 - Extended Rear legs required p/n 020-10016-00 for models 6211, 6214, 6215-001, 6215-002, 6215-003, 6217



Midmark® 6215-002

Compact Rectangular Workstation

- 26" of height adjustment from 23" to 49"
- 17.5" D x 23.5" W worksurface
- 19.5" D x 20.5" W base dimension
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Weight capacity: 20 lb
- Unit weight: 48 lb
- Shipping weight: 93 lb

Model

6215-002-xxx¹

Description

Compact Workstation

Accessories

Plant Installed

Field Installed

BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-00 ²	Equipment Pole Mount for Digital Vital Signs Device
MPSF	9A648001	Mobile Printer Shelf
PMSL ³	9A655001 ³	Power Mount w/ Steel Cable Mgmt-LH
PMSR ³	9A655001 ³	Power Mount w/ Steel Cable Mgmt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Power Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Power Cord-RH
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB ⁴	9A656001 ⁴	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount

¹ - See color code chart on page 24

² - Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.

³ 9A655001 & 9A655002 - Extended Rear legs required p/n 020-10016-00 for models 6211, 6214, 6215-001, 6215-002, 6217

⁴ 9A656001 - Extended Rear legs required p/n 020-10016-00 for models 6211, 6214, 6215-001, 6215-002, 6215-003, 6217

MTR CODE

HAW4



Midmark® 6215-003

Compact Kidney Workstation

- 26" of height adjustment from 23" to 49"
- 17" D x 25" W worksurface
- 19.5" D x 20.5" W base dimension
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Weight capacity: 20 lb
- Unit weight: 48 lb
- Shipping weight: 93 lb

Model

6215-003-xxx¹

Description

Compact Workstation

Accessories

Plant Installed

BCHL
BCHR
CPUB
CPUF
DBCM
SBCB
SBCF
SCML
SCMR
SHLH
SHRH
USCB²

Field Installed

9A659001
9A659001
9A663001
9A663001
9A653001
9A652001
9A652001
9A668001
9A668001
9A654001
9A654001
9A656001²

Bar Code Holder-LH
Bar Code Holder-RH
CPU Holder-Back
CPU Holder-Front
Double Basket
Single Basket-Back
Single Basket-Front
Steel Brkt + MG Coiled Power Cord-LH
Steel Brkt + MG Coiled Power Cord-RH
Scale Holder-LH
Scale Holder-RH
Adjustable Utility Shelf

¹ - See color code chart on page 24

Requires 9A651001 or 9A658001 for field installed option.

² 9A656001 - Extended Rear legs required p/n 020-10016-00 for models 6211, 6214, 6215-001, 6215-002, 6215-003, 6217



Midmark® 6217

Secure Laptop Workstation + Mouse Tray

- 26" of height adjustment from 23" to 49"
- Keyed laptop security system
- 23" D x 22" W worksurface w/ slideout mouse surface
- 11" D x 20.5" W secondary worksurface
- 19.5" D x 20.5" W base dimension
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Weight capacity: 30 lb
- Unit weight: 60 lb
- Shipping weight: 105 lb

Model

6217-001-xxx¹

Description

Secure Laptop Workstation + Mouse Tray

Accessories

Plant Installed

Field Installed

BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/ Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-00 ²	Equipment Pole Mount for Digital Vital Signs Device
MPSF	9A648001	Mobile Printer Shelf
PMSL ³	9A655001 ³	Power Mount w/ Steel Cable Mgmt-LH
PMSR ³	9A655001 ³	Power Mount w/ Steel Cable Mgmt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Power Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Power Cord-RH
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB ⁴	9A656001 ⁴	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount

¹ See color code chart on page 24

² Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.

³ 9A655001 & 9A655002 - Extended Rear legs required p/n 020-10016-00 for models 6211, 6214, 6215-001, 6215-002, 6217

⁴ 9A656001 - Extended Rear legs required p/n 020-10016-00 for models 6211, 6214, 6215-001, 6215-002, 6215-003, 6217



Midmark® 6218

Flat Panel PC Workstation + Mouse Tray

- 26" of height adjustment from 23" to 49"
- Monitor arm capacity range 5 lb min - 20 lb max
- Monitor arm independent rotation of 180°
- 23" D x 22" W worksurface w/ slideout mouse surface
- 22.5" D x 20.5" W base dimension
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical adj. and up to 40° positive tilt
- Weight capacity: 30 lb
- Unit weight: 55 lb
- Shipping weight: 100 lb

Model

6218-001-xxx¹

Description

Flat Panel PC Workstation + Mouse Tray

Accessories

Plant Installed

Field Installed

AIOM	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLM	9A647001	Drawer w/ Lock-Mouse Trk Mnt
DWNM	9A646001	Drawer-Mouse Track Mnt
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-00 ²	Equipment Pole Mount for Digital Vital Signs Device
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
PMSL ³	9A655001 ³	Power Mount w/ Steel Cable Mgmt-LH
PMSR ³	9A655001 ³	Power Mount w/ Steel Cable Mgmt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Power Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Power Cord-RH
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB ⁴	9A656001 ⁴	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount

¹ See color code chart on page 24

² Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.



Midmark® 6219

Flat Panel Secure PC Workstation + Mouse Tray

- 26" of height adjustment from 23" to 49"
- Monitor arm capacity range 5 lb min - 20 lb max
- Keyed PC security system
- Monitor arm independent rotation of 180°
- 23" D x 22" W worksurface w/ slideout mouse surface
- 22.5" D x 20.5" W base dimension
- 11" D x 20.5" W secondary worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical adj. + up to 40° positive tilt
- Weight capacity: 30 lb
- Unit weight: 64 lb
- Shipping weight: 109 lb

Model

6219-001-xxx¹

Description

Flat Panel Secure PC Workstation + Mouse Tray

Accessories

Plant Installed

Field Installed

AIOM	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLM	9A647001	Drawer w/ Lock-Mouse Trk Mnt
DWNM	9A646001	Drawer-Mouse Track Mnt
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-00 ²	Equipment Pole Mount for Digital Vital Signs Device
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
PCML	9A669001	MG Coiled Power Cord + Steel Hanger - LH
PCMR	9A669001	MG Coiled Power Cord + Steel Hanger - RH
PMBL ³	9A655002 ³	Power Mount w/ Bin Cable Mgmnt-LH
PMBR ³	9A655002 ³	Power Mount w/ Bin Cable Mgmnt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB ⁴	9A656001 ⁴	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount

¹ See color code chart on page 24

² Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.



Midmark® 6281

Workstation

- 22" height adjustment
- 19.5" W x 14" D worksurface
- Slide-out mouse surface
- Up to 45" reach off wall
- Weight capacity: 0-8 lb

- Folds to 9" flat to wall
- Unit weight: 33 lb
- Shipping weight: 40 lb
- Work surface accommodates up to 17" laptops

MTR CODE
HAW4 CS4

Model	Description
6281-001	Laptop Wall Mounted Workstation

Accessories

Plant Installed	Field Installed	
n/a	9A588001	20" Extension Arm
n/a	9A591001	Track Mount Bar Code Holder
n/a	9A590001	Tamper Proof PC Holder
n/a	9A589001	Adjustable PC Holder

MTR CODE

HAW4
CS4



Midmark® 6282 Workstation

- 22" worksurface height adjustment
- 25" W x 10" D worksurface
- Max monitor capacity is 20 lb
- 9.5" independent monitor height adjustment
- Independent 180° monitor rotation and 0-45° tilt
- VESA® mount 75 mm and 100 mm
- Folds to 9.47" flat to wall
- Weight capacity: 5 - 20 lb
- Product weight: 43 lb
- Shipping weight: 50 lb

Model

6282-001

Description

Standard Duty Wall Mounted Workstation

Accessories

Plant Installed

n/a
n/a
n/a
n/a
n/a
n/a
n/a

Field Installed

9A593001
9A588001
9A592001
9A591001
9A590001
9A589001
9A660001

Dual Monitor Arm
20" Extension Arm
VESA Mount Bar Code Holder
Track Mount Bar Code Holder
Tamper Proof PC Holder
Adjustable PC Holder
VESA Mounting Plate

MTR CODE

HAW4
CS4



Midmark® 6283

Workstation

- 22" worksurface height adjustment
- 25" W x 10" D worksurface
- Max monitor capacity is 35 lb
- 9.5" independent monitor height adjustment
- Independent 180° monitor rotation and 0-45° tilt
- VESA® mount 75 mm and 100 mm
- Folds to 9.47" flat to wall
- Weight capacity: 20-35 lb
- Product weight: 43 lb
- Shipping weight: 50 lb

Model

6283-001

Description

Heavy Duty Wall Mounted Workstation

Accessories







Plant Installed	Field Installed	
n/a	9A593001	Dual Monitor Arm
n/a	9A588001	20" Extension Arm
n/a	9A592001	VESA Mount Bar Code Holder
n/a	9A591001	Track Mount Bar Code Holder
n/a	9A590001	Tamper Proof PC Holder
n/a	9A589001	Adjustable PC Holder
n/a	9A660001	VESA Mounting Plate

Midmark Synthesis® Mobile Solutions are designed for the clinical space to be flexible and durable to stand the test of time.






Frame/Panel Colors:














Solids:

	216	Pebble Grey
	793	Pebble Grey PVC-Free
	490	Fawn
	649	Whisper
	791	Frost
	931	Frost PVC Free

Metallics:

	804	Venus Silver (M)
	805	Pearl Essence (M)
	806	Earthen Bronze (M)
	807	Flat Iron (M)

Woodgrains:

	459	Acorn
	487	Earth
	794	Dune
	795	Flax
	796	Nest
	797	Hewn
	798	Path
	799	Sediment
	800	Radiance
	801	Timber
	802	Storm
	803	Henna
	874	Brazilwood



* Shown with Kydex®
contour top

Mobile Treatment Cabinet M2

33" H x 27" W x 18" D (2) 4" Drawers + (2) Doors

- Serenity panel style w/ flare handles
- Shipping Weight: 155 lb

Product Information

Step 1: Configuration

M2

Step 2: Depth

1 18" Deep

Step 3: Option Package

- A 60 mm Casters
- B 4" Quiet Casters
- C 60 mm Casters + Keyed Locks
- D 4" Quiet Casters + Keyed Locks

Step 4: Cabinetry Finish*

- XXX Solid Color Cabinet Finish (Excludes PVC-Free)
- XXX Woodgrain, Metallic + PVC-Free Cabinet Finish

Step 5: Countertop Type

- KCN Kydex Contour Top
- SEB Soft Edge Bumper Top
- SEH Soft Edge Handle Top (29"W)

Configuration Example: M21A.791.KCN

*See page 24 for Cabinetry Finish

Kydex is a registered trademark of SEKISUI POLYMER INNOVATIONS, LLC



Mobile Treatment Cabinet M4

33" H x 27" W x 18" D (3) 4" Drawers + (1) 8" Drawer

- Serenity panel style w/ flare handles
- Shipping Weight: 152 lb

* Shown with soft-edge bumper top

Product Information

Step 1: Configuration

M4

Step 2. Depth

1 18" Deep

Step 3. Option Package

- A** 60 mm Casters
- B** 4" Quiet Casters
- C** 60 mm Casters + Keyed Locks
- D** 4" Quiet Casters + Keyed Locks

Step 4. Cabinetry Finish*

- XXX** Solid Color Cabinet Finish (Excludes PVC-Free)
- XXX** Woodgrain, Metallic + PVC-Free Cabinet Finish

Step 5. Countertop Type

- KCN** Kydex® Contour Top
- SEB** Soft Edge Bumper Top
- SEH** Soft Edge Handle Top (29"W)

Configuration Example: M41A.791.SEB

*See page 24 for Cabinetry Finish

Kydex is a registered trademark of SEKISUI POLYMER INNOVATIONS, LLC



Mobile Treatment Cabinet M5

33" H x 27" W x 18" D (5) 4" Drawers

- Serenity panel style w/ flare handles
- Shipping Weight: 159 lb

* Shown with soft-edge handle top

Product Information

Step 1: Configuration

M5

Step 2. Depth

1 18" Deep

Step 3. Option Package

- A** 60 mm Casters
- B** 4" Quiet Casters
- C** 60 mm Casters + Keyed Locks
- D** 4" Quiet Casters + Keyed Locks

Step 4. Cabinetry Finish*

XXX Solid Color Cabinet Finish (Excludes PVC-Free)

XXX Woodgrain, Metallic + PVC-Free Cabinet Finish

Step 5. Countertop Type

- KCN** Kydex® Contour Top
- SEB** Soft Edge Bumper Top
- SEH** Soft Edge Handle Top (29"W)

Configuration Example: M51A.791.SEH

*See page 24 for Cabinetry Finish

Kydex is a registered trademark of SEKISUI POLYMER INNOVATIONS, LLC



* Shown in 24" Width

Mobile Procedure Cart P1

- 43" H x 18" D and available in 18", 21" and 24" W
- (2) 2" drawers, (3) 4" drawers and (1) 8" drawer
- Standard features: slide-out worksurface, no-tip base, full extension Drawers, 4" casters and base bumper
- Various accessories and accessory mounting locations
- Cascade panel style w/ flare handles
- Shipping Weight:
 - o 18" – 187 lb
 - o 21" – 199 lb
 - o 24" – 205 lb

Product Information

Step 1: Configuration

P1 (2) 2" drawers, (3) 4" drawers, (1) 8" drawer
Slide-out Worksurface, 4" Casters, Base Bumper

Step 2. Width

- 1** 18" wide
2 21" wide
3 24" wide

Step 3. Options

- N** None
E Keyed Lock
F Keyless Lock

Step 4. Cabinetry Finish

- XXX** Solid Color Cabinet Finish (Excludes PVC-Free)
XXX Woodgrain, Metallic + PVC-Free Cabinet Finish

Step 5. Countertop Type

- KCT** Kydex® Contour Top
SCT Soft-Edge Contour Top

Step 6. Accessories

Plant Installed Field Installed LHHE, RHHE, LHUA, LHDA, RHUA, RHDA and RPCK are PLANT INSTALLED ONLY

LHHE	---	Side Transit Handle, LH mount
RHHE	---	Side Transit Handle, RH mount
LHUA	---	Side Accessory Rail, LH mount - upper only
LHDA	---	Side Accessory Rails, LH mount - upper and lower rails
RHUA	---	Side Accessory Rail, RH mount - upper only
RHDA	---	Side Accessory Rails, RH mount - upper and lower rails
RBCK	---	Top Rear Accessory Rack (21" + 24" carts hold 3 accessories, 18" holds 2)

Step 7. These accessories below require one or more of LHUA, LHDA, RHUA, RHDA or RBCK to hold them

Plant Installed	Field Installed	
SBL	9A688001	6 x 6 Storage Storage Bins/Hinged Lid
MWC	9A689001	Waste Basket w/Bracket (fits lower side accessory rail only)
SWH	9A690001	Disinfectant Wipe Bracket
MSC	9A691001	Sharps container w/Bracket (fits upper side accessory rail only)
SGH	9A692001	Single Glove Box Holder
TGH	9A693001	3-Box Glove Holder (fits top or side on 21" + 24" cart, side only on 18")

Configuration Example: P12E.216.KCT (Does not include accessories)



* Shown in 18" Width

Mobile Procedure Cart P2

- 43" H x 18" D and available in 18", 21" and 24" W
- (4) 4" drawers and (1) 8" drawer
- Standard features: slide-out worksurface, no-tip base, full extension Drawers, 4" casters and base bumper
- Various accessories and accessory mounting locations
- Cascade panel style w/ flare handles
- Shipping Weight:
 - o 18" – 187 lb
 - o 21" – 199 lb
 - o 24" – 205 lb

Product Information

Step 1: Configuration

P2 (4) 4" drawers, (1) 8" drawer
Slide-out Worksurface, 4" Casters, Base Bumper

Step 2: Width

1 18" wide
2 21" wide
3 24" wide

Step 3: Options

N None
E Keyed Lock
F Keyless Lock

Step 4: Cabinetry Finish

XXX Solid Color Cabinet Finish (Excludes PVC-Free)
XXX Woodgrain, Metallic + PVC-Free Cabinet Finish

Step 5: Countertop Type

KCT Kydex® Contour Top
SCT Soft-Edge Contour Top

Step 6: Accessories

Plant Installed Field Installed LHHE, RHHE, LHUA, LHDA, RHUA, RHDA and RPKC are PLANT INSTALLED ONLY

LHHE	---	Side Transit Handle, LH mount
RHHE	---	Side Transit Handle, RH mount
LHUA	---	Side Accessory Rail, LH mount - upper only
LHDA	---	Side Accessory Rails, LH mount - upper and lower rails
RHUA	---	Side Accessory Rail, RH mount - upper only
RHDA	---	Side Accessory Rails, RH mount - upper and lower rails
RBCK	---	Top Rear Accessory Rack (21" + 24" carts hold 3 accessories, 18" holds 2)

Step 7: These accessories below require one or more of LHUA, LHDA, RHUA, RHDA or RBCK to hold them

Plant Installed	Field Installed	
SBL	9A688001	6 x 6 Storage Storage Bins/Hinged Lid
MWC	9A689001	Waste Basket w/Bracket (fits lower side accessory rail only)
SWH	9A690001	Disinfectant Wipe Bracket
MSC	9A691001	Sharps container w/Bracket (fits upper side accessory rail only)
SGH	9A692001	Single Glove Box Holder
TGH	9A693001	3-Box Glove Holder (fits top or side on 21" + 24" cart, side only on 18")

Configuration Example: P23F.794.SCT (Does not include accessories)

*See page 24 for Cabinetry Finish

Kydex is a registered trademark of SEKISUI POLYMER INNOVATIONS, LLC



* Shown in 21" Width

Mobile Procedure Cart P3

- 43" H x 18" D and available in 18", 21" and 24" W
- (2) 4" drawers, (1) 23" door w/1 adjustable shelf
- Standard features: slide-out worksurface, no-tip base, full extension Drawers, 4" casters and base bumper
- Various accessories and accessory mounting locations
- Cascade panel style w/ flare handles
- Shipping Weight:
 - o 18" – 187 lb
 - o 21" – 199 lb
 - o 24" – 205 lb

Product Information

Step 1: Configuration

P3 (2) 4" drawers, (1) 23" door w/1 adjustable shelf
Slide-out Worksurface, 4" Casters, Base Bumper

Step 2: Width

- 1** 18" wide (RH single door)
2 21" wide (RH single door)
3 24" wide (double doors)

Step 3: Options

- N** None
E Keyed Lock
F Keyless Lock

Step 4: Cabinetry Finish

- XXX** Solid Color Cabinet Finish (Excludes PVC-Free)
XXX Woodgrain, Metallic + PVC-Free Cabinet Finish

Step 5: Countertop Type

- KCT** Kydex® Contour Top
SCT Soft-Edge Contour Top

Step 6: Accessories

Plant Installed Field Installed LHHE, RHHE, LHUA, LHDA, RHUA, RHDA and RPCK are PLANT INSTALLED ONLY

LHHE	---	Side Transit Handle, LH mount
RHHE	---	Side Transit Handle, RH mount
LHUA	---	Side Accessory Rail, LH mount - upper only
LHDA	---	Side Accessory Rails, LH mount - upper and lower rails
RHUA	---	Side Accessory Rail, RH mount - upper only
RHDA	---	Side Accessory Rails, RH mount - upper and lower rails
RBCK	---	Top Rear Accessory Rack (21" + 24" carts hold 3 accessories, 18" holds 2)

Step 7. These accessories below require one or more of LHUA, LHDA, RHUA, RHDA or RBCK to hold them

Plant Installed	Field Installed	
SBL	9A688001	6 x 6 Storage Storage Bins/Hinged Lid
MWC	9A689001	Waste Basket w/Bracket (fits lower side accessory rail only)
SWH	9A690001	Disinfectant Wipe Bracket
MSC	9A691001	Sharps container w/Bracket (fits upper side accessory rail only)
SGH	9A692001	Single Glove Box Holder
TGH	9A693001	3-Box Glove Holder (fits top or side on 21" + 24" cart, side only on 18")

Configuration Example: P31E.791.KCT (Does not include accessories)

MSS1
MSS2
PP1
SC1
SH1



* 6" x 12" Storage Bin Shown

Midmark® Storage Bins

- Clinical space organization and storage
- Works great with Midmark angled flow-shelving
- Stackable bins, translucent lids, bin dividers

Model	Description
002-10708-00	6" x 6" Storage Bins, 8 ct box
002-10709-00	3" x 12" Storage Bins, 4 ct box
002-10710-00	6" x 12" Storage Bins, 4 ct box
002-10711-00	12" x 12" Storage Bins, 2 ct box
002-10712-00	6" x 6" Storage Bin Lids, 8 ct pack
002-10713-00	3" x 12" Storage Bin Lids, 4 ct pack
002-10714-00	6" x 12" Storage Bin Lids, 4 ct pack
002-10715-00	12" x 12" Storage Bin Lids, 2 ct pack
002-10716-00	6" Storage Bin Dividers, 4 ct pack
002-10717-00	12" Storage Bin Dividers, 4 ct pack



Ritter® Waste Receptacle

- Includes bag retaining ring
- 14" W x 16" D x 25" H
- Shipping weight: 17 lb
- Capacity: 43 liquid qts
- Color: Pebble

Model

261-001

Description

Biohazard waste receptacle

Glove, Towel + Cup Dispenser

- D5 – D6 = 18"W x 9"D x 24"H (Vertical)
- D7 – D8 = 24"W x 9"D x 24"H (Vertical)
- Shipping weight: 37 lb
- Serenity panel style w/flare handle



* D6 shown
Vertical



* D8 shown
Vertical

Product Information

Step 1. Configuration

- D5** 18" Vertical Glove/Towel/Cup Dispenser - (2) Opening RH Door
- D6** 18" Vertical Glove/Towel/Cup Dispenser - (2) Opening LH Door
- D7** 24" Vertical Glove/Towel/Cup Dispenser - (3) Opening RH Door
- D8** 24" Vertical Glove/Towel/Cup Dispenser - (3) Opening LH Door

Step 2. Option Package

- X** No Lock

Step 3. Cabinetry Finish*

- XXX** Solid Color Cabinet Finish (Excludes PVC-Free)
- XXX** Woodgrain, Metallic + PVC-Free Cabinet Finish

Configuration Example: D5X.803

*See page 24 for Cabinetry Finish



* N2 shown

Dressing Nook

- Right-hand and left-hand door pictured
- 24" W x 7" D x 63" H
- Interior mirror, garment hooks and shoe storage included
- Serenity panel style w/flare handle
- Shipping weight: 138 lb



* N1 shown

Product Information

Step 1. Configuration

N1 Dressing Nook - RH Door

N2 Dressing Nook - LH Door

Step 2. Option Package

X No Lock

A Keyed Locks


Step 3. CabinetryFinish*

XXX Solid Color Cabinet Finish (Excludes PVC-Free)

XXX Woodgrain, Metallic + PVC-Free Cabinet Finish

Configuration Example: N2X.803

*See page 24 for Cabinetry Finish



Click here to visit:
Synthesis® Cabinetry Literature

MTR CODE
MSS1
MSS2
PP1
SC1
SH1



Click here to visit:
Synthesis® Cabinetry Style Guide

MTR CODE
MSS1
MSS2
PP1
SC1
SH1

Midmark® + Ritter® Medical Equipment (Tables, Chairs, Physician Seating, Lights and Sterilizers) Limited Warranty

SCOPE OF WARRANTY: Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will, at Midmark's option, repair or replace components of the domestic and international medical products manufactured by Midmark (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. The sole remedy under this limited warranty is the repair or replacement, at Midmark's option, of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which are determined to exist upon examination by Midmark. This warranty extends only to the original retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

Midmark warrants to the original retail purchaser that during the applicable warranty period it will repair or replace software contained within the products manufactured by Midmark (except for those not warranted under "Exclusions") if: (1) the media on which the software is furnished exhibits defects in material or workmanship under normal use; or (2) the software does not substantially conform to its published specifications.

APPLICABLE WARRANTY PERIOD: The applicable warranty period, measured from the date of the original purchase of the retail product shall be:

- One (1) year for all standard warranted products and components
- Two (2) years for Midmark M9/M11-05X Steam Sterilizer models
- Three (3) years for QuickClean® Ultrasonic Cleaners

OBTAINING WARRANTY SERVICE: Warranty service must be obtained through either Midmark or an authorized dealer in the Midmark product line for which warranty service is requested. Midmark may be contacted for warranty service inquiries or issues via email at midmark.com; by phone at 1.800.MIDMARK or by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380.

It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this limited warranty.

EXCLUSIONS: This limited warranty does not cover and Midmark shall not be liable for the following:

- (1) Defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) Products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environmental conditions and electrical requirements;
- (3) Products considered to be of a consumable or sterile nature;
- (4) Accessories or parts not manufactured by Midmark;
- (5) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (6) Costs and expenses of routine maintenance and cleaning; and
- (7) Representations and warranties made by any person or entity other than Midmark.
- (8) Matching of color, grain or texture except to commercially acceptable standards;
- (9) Changes in color caused by natural or artificial light;
- (10) Custom manufactured products;
- (11) Alterations or modifications to the product by any person or entity other than Midmark; and
- (12) Products that would otherwise be covered under this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical market.

SOFTWARE; WITH RESPECT TO SOFTWARE THAT IS A PRODUCT OR COMPONENT THEREOF, MIDMARK DOES NOT WARRANT THAT THE SOFTWARE: (1) IS ERROR FREE; (2) CAN BE USED WITHOUT PROBLEMS OR INTERRUPTIONS; OR (3) IS FREE FROM VULNERABILITY TO INTRUSION OR ATTACK BY VIRUSES OR OTHER METHODS.

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER: MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, LOSS OF DATA, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. THIS DISCLAIMER SHALL SURVIVE ANY FAILURE OR ASSERTED FAILURE OF THE ESSENTIAL PURPOSE OF THIS LIMITED WARRANTY OR ITS REMEDIES SPECIFIED HEREIN.

WARRANTY DISCLAIMER: THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

STATUTE OF LIMITATIONS: No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

NO AUTHORIZATION: No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

Warranty information included in product manuals.

Midmark® Workstations: Limited Warranty

SCOPE OF WARRANTY

Midmark Corporation warrants to the original purchaser its Midmark Workstation products manufactured by Midmark to be free from defects in material and workmanship under normal interior use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of the parts or the products the defects of which are reported to Midmark within the applicable warranty period and which upon examination by Midmark prove to be defective.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows:

- A. Five years for all Midmark Workstation products and components (except for the shorter periods in B), including drawer fronts and casters
- B. One year for power system components and coiled power cords

EXCLUSIONS

This warranty does not cover and Midmark shall not be liable for: (1) Failure to install or use according to the Midmark user and installation manuals (2) Normal wear and tear, repairs and replacements required because of misuse, abuse, negligence, alteration, accident, freight damage or tampering (including, without limitation, use of product in unsuitable environments or conditions); (3) matching of color, grain or texture except to commercially acceptable standards; (4) changes in color caused by natural or artificial light; (5) products which are

not installed, used or properly cleaned as required in Midmark's written Installation Manual; (6) specially manufactured products; (7) products considered to be of a consumable nature such as light bulbs and surge suppression products; (8) accessories or parts not manufactured by Midmark or use of non-authorized Midmark components; (9) charges by anyone (including Midmark's authorized dealers) for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which is not expressly authorized in writing in advance by Midmark.

EXCLUSIVE REMEDY

Midmark's only obligation under this warranty is the repair or replacement of defective parts. Midmark shall not be liable for any direct, special, indirect, incidental, exemplary or consequential damages or delay.

NO AUTHORIZATION

No person or firm is authorized to create for Midmark any other obligation or liability in connection with the products.

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

Synthesis® Cabinetry Collection by Midmark: Limited Warranty

SCOPE OF WARRANTY

Midmark Corporation warrants to the original purchaser its Midmark Modular cabinetry products (Synthesis cabinetry Collection) manufactured by Midmark to be free from defects in material and workmanship under normal interior use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of the parts or the products the defects of which are reported to Midmark within the applicable warranty period and which upon examination by Midmark prove to be defective.

Note: If the Midmark Modular Cabinetry products are installed by Midmark, the above referenced warranty shall also include any labor performed, which has been authorized in writing in advance by Midmark, in connection with such repair or replacement during the applicable warranty period.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows:

- A. Ten years for all Midmark brand cabinetry products and components (except for the shorter periods in B and C), including door and drawer fronts and casters*
- B. Three years for task lights
- C. One year for countertop surfaces and accessories

*The new 10-year warranty will apply to original users who took delivery of the Midmark brand cabinetry products and components in question on or after January 1, 2023.

SERVICE PERIOD

Ten years, with exceptions on discontinuation of colors, components and parts.

If a color is discontinued, Midmark will service it for the longer of 5-years or until supplies are no longer available, at which time a current color must be selected.

On components and parts: in the event of a discontinuation, Midmark will offer a comparable replacement component or part.

EXCLUSIONS

This warranty does not cover sinks, faucets and plumbing accessories. Any warranties on these items are extended directly by the manufacturer of these items to the original purchaser. Information on these manufacturers' warranties will be enclosed with the applicable products. Also, Midmark will furnish copies of any warranties extended by any such manufacturers upon request.

This warranty does not cover and Midmark shall not be liable for: (1) repairs and replacements required because of misuse, abuse, negligence, alteration, accident, freight damage or tampering; (2) matching of color, grain or texture except to commercially acceptable standards; (3) changes in color caused by natural or artificial light; (4) products which are not installed, used or properly cleaned as required in Midmark's written Installation Manual; (5) specially manufactured products; (6) products considered to be of a consumable nature such as light bulbs and surge suppression products; (7) accessories or parts not manufactured by Midmark; (8) charges by anyone (including Midmark's authorized dealers) for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which is not expressly authorized in writing in advance by Midmark.

EXCLUSIVE REMEDY

Midmark's only obligation under this warranty is the repair or replacement of defective parts. Midmark shall not be liable for any direct, special, indirect, incidental, exemplary or consequential damages or delay.

NO AUTHORIZATION

No person or firm is authorized to create for Midmark any other obligation or liability in connection with the products.

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.



Designing better care.®



CARB 93120.2 Phase 2 Compliant and TSCA Title VI Compliant
(applicable products)

The color examples shown are the best representation of the original material. Actual color may vary slightly. We strongly recommend that you contact Midmark Customer Experience at 1.800.MIDMARK to request a sample before placing your order.

Midmark is an ISO 13485 and ISO 9001 Certified Company.
Certain products are not included. See the complete list at: midmark.com/ISO

For more information, contact your Midmark dealer
or call: 1.800.MIDMARK Outside the USA call: 1.937.526.3662
or visit our website: midmark.com

© 2025 Midmark Corporation, Versailles, Ohio USA
Products subject to improvement changes without notice.

007-10523-00 (2/25)

