

# VetPro® 1000 Antimicrobial Cover 9A477001

## ***Recommended Easy Cleaning..***

- A) To remove soils or stains, wipe clean with soapsuds and lukewarm water.*
- B) To remove hard to clean spots, use a brush with a standard liquid household cleaner. If necessary, pre-soak. Do not use solvents or harsh cleaners.*
- C) Phenolic or quaternary type disinfectants can be used with cover when diluted per manufacturer's instructions. Iodophor type disinfectants are not recommended for use on cover.*

Midmark Corporation  
60 Vista Drive  
Versailles, OH 45380 USA  
1-800-643-6275  
1-937-526-3662

