



MA353102

**9A03004 Armboard Assembly w/Extra Bracket**  
(includes items 1 thru 11)

**9A03005 Drawer Spacer** (incl. items 12 & 13)  
(required for tables w/drawers)

Item	Description	Qty.
1	Armrest Board	1
2	Washer	2
3	Bracket	1
4	Washer	1
5	Cap Nut	1
6	Hex Nut	1
7	Knob	2
8	Bracket	2
9	Screw	4
10	Screw	4
11	Plate Assembly	1
12	Spacer Weldment (for 203 models w/drawers)	2
13	Screw (#10-24 x 1 1/4")	6

Serial Numbers / Dates: all