

## Specifications / Range of Motion: 6251 / 6252 / 6256

Specifications: 6251 / 6252 / 6256	
<b>Lift Capacity:</b>	Lift capacity varies depending on the lift system, and accessories that are installed.  Carts are configurable with <i>Standard Duty</i> , or <i>Heavy Duty</i> lift systems.  Standard Duty: 13-40 lbs. Heavy Duty: 41-60 lbs.
<b>Base Dimensions:</b>	18" x 20" (46 cm x 51 cm)
<b>Electrical System Information:</b>	
Input Power Requirements	Rated Voltage: AC 100-240V, 50/60 Hz Rated Current: 6.4 A max.
Output Power	DC 10.5v-15V (Load) 20A max, 300W max 3.0-15V (battery)
<b>Battery Information:</b>	
Sealed Lead Acid Battery (SLA)	<b>Size:</b> 55 Amp/Hr. (AH) <b>Charging:</b> up to 5 hours for 100% <b>Cycle Life:</b> 200 to 1,000 recharge cycles (approx. 18 to 36 months)
<b>Worksurface Dimensions:</b>	
<b>6251 / 6252</b>	20" D x 20" W (51 cm x 51 cm)
<b>6256</b>	22" D x 32" W (56 cm x 81 cm)
<b>Regulatory Compliance:</b>	ANSI / BIFMA UL 60601-1 CAN/CSA-C22.2 No. 601-1-M90

