

Specifications / Range of Motion: 6261 / 6263

Specifications: 6261 / 6263	
Lift Capacity:	Lift capacity varies depending on the lift system, and accessories that are installed. Carts are configurable with <i>Standard Duty</i> , or <i>Heavy Duty</i> lift systems. Standard Duty: 13-40 lbs. Heavy Duty: 41-60 lbs.
Base Dimensions:	18" x 20" (46 cm x 51 cm)
Electrical System Information:	
Input Power Requirements	120 V, 60 Hz, 5.0A
Output Power	120 V, 60 Hz 150 VA, 150 W (nominal)
Battery Information:	
Lithium Iron Phosphate Battery (LiFe)	Size: 36.8 Amp/Hr. (AH) Charging: 100% recharged in 2 hours Cycle Life: 8,000 cycles / 5 years Weight: 11 lbs. (5 kg)
Sealed Lead Acid Battery (SLA)	Size: 55 Amp/Hr. (AH) Charging: up to 5 hours for 100% Cycle Life: up to 1,200 cycles Weight: 39 lbs. (17.7 kg)
Worksurface Dimensions:	Lower Worksurface: 20" W x 20" D (51 x 51 cm) Upper Worksurface: 20" W x 9" D (51 x 23 cm)
Regulatory Compliance:	ANSI / BIFMA IEC 60601-1 UL 60601 UL 1778

