




Notes:



Leviton 42080-1xS or equivalent
Leviton 41106-Rx6 or equivalent RJ11



Leviton 42080-2xS or equivalent
Leviton 41106-Rx6 or equivalent RJ11
Leviton 5G108-Rx5 or equivalent RJ45



Leviton 42079-1xS or equivalent
Leviton 5G108-Rx5 or equivalent RJ45

* Total length of wiring between devices, including cabling outside the wall.

**Midmark
IQvitals Zone®/626
Pre-wire Diagram
004-10045-00**