




Service Bulletin - RoHS Compliant PC Boards

Note

RoHS (Restriction of Hazardous Substances) compliant PC Boards are made of more environmentally responsible materials and still offer the same functionality as our existing PC boards.

Midmark will be using RoHS compliant PC Boards in our new products released in 2014. The new boards can be easily identified by their RED color.

RoHS compliant PC Boards will also be gradually used for service and replacement parts as well. Part numbers of the existing PC boards will not change and can be used interchangeably with all existing green colored boards. Products having multiple PC boards can use a mix of the newer red with the existing green boards without experiencing any problems.



Note

Red boards will replace green boards. Their functionality is equal to and fully interchangeable with existing PC Boards.

TP201 Rev. A