
























Midmark Wall Mounted Workstation, Service Parts List

Dec 17-2015

Part Number	Description	Image	Laptop	HD	LT	Number Needed per Assembly	
1210268	BASE SWIVEL BUSHING SLEEVE		Yes	Yes	Yes	1	
1210269	ARM FRONT SWIVEL BUSHING		Yes	Yes	Yes	1	
1210274	END OF ARM BOLT 3/8-16 X 5.75 BHCS LOCK PATCH BLK		Yes	Yes	Yes	1	
1210371	MONITOR TOWER SPRING ATTACHMENT PINS		Yes	Yes	Yes	3	
1210577	MAIN ARM CABLE MANAGER SCREWS		Yes	Yes	Yes	4	
1330601	WALL MOUNTING KIT (5 X SCREWS 1/4-20, 10 X SCREWS M6, 10 X TOGGLERS)		Yes	Yes	Yes	1	
1330602	ARM RETENTION HARDWARE KIT (BOLT, WASHERS, PINS)		Yes	Yes	Yes	1	
6011616	VESA QUICK CONNECT BRACKET, MONITOR MOUNTING SCREWS, CABLE MANAGER CLIP AND VELCRO KIT		No	Yes	Yes	1	
6021150	DUAL MONITOR RAIL ADJUSTMENT KNOB		No	Yes	Yes	2	
7000168	MONITOR TOWER TOP CAP SCREW (#10 X .250 BHCS)		No	Yes	Yes	2	
7902173	WORKSTATION MOUNTING SCREW		No	Yes	Yes	2	
120045162	WRIST REST (FOR MIDMARK WORKSTATION)		No	Yes	Yes	1	
1210267RX	MAIN ARM SIDE SHROUD, RIGHT HAND		Yes	Yes	Yes	1	






Midmark Wall Mounted Workstation, Service Parts List

Dec 17-2015

Part Number	Description	Image	Laptop	HD	LT	Number Needed per Assembly	
1100164W	SPIRAL WRAP FOR CABLES, WHITE		Yes	Yes	Yes	1	
1210260LX	MONITOR TOWER TILT SHROUD , LEFT HAND, GREY		No	Yes	Yes	1	
1210260RX	MONITOR TOWER TILT SHROUD , RIGHT HAND, GREY		No	Yes	Yes	1	
1210262AG	MAIN ARM ADJUSTMENT COVER CAP, GREY		Yes	Yes	Yes	1	
1210265X	MAIN ARM CABLE MANAGER, GREY		Yes	Yes	Yes	1	
1210267LX	MAIN ARM SIDE SHROUD, LEFT HAND		Yes	Yes	Yes	1	
1210361X	MONITOR TOWER TOP CAP, GREY		No	Yes	Yes	1	
1210364X2400	CABLE COVER FOR TRACK, GREY		Yes	Yes	Yes	1	
1210368X	KEYBOARD TILT SIDE CAPS, GREY		Yes	Yes	Yes	1	
1210600X	MAIN ARM ASSEMBLY, COMPLETE, GREY		Yes	Yes	Yes	1	
1210602E	MONITOR TILT SUBASSEMBLY - BLACK SATIN		No	Yes	Yes	1	
1210607X	MONITOR TOWER, COMPLETE, FREY		No	Yes	Yes	1	
1210611X	WALL EAR ASSEMBLY, GREY		Yes	Yes	Yes	1	

Midmark Wall Mounted Workstation, Service Parts List

Dec 17-2015

Part Number	Description	Image	Laptop	HD	LT	Number Needed per Assembly	
130050162X	TOP CAP FOR LAPTOP-ARM JOINT, GREY		Yes	No	No	1	
133050170	LAPTOP MOUNTING SCREWS		Yes	No	No	4	
133050MD600	LAPTOP HARDWARE SERVICE KIT (BOLT, WASHERS, BORE CUP, TOP CAP) PANTONE PQ COOL GRAY 1C ("MIST")		Yes	No	No	1	
1330603X	EXTENSION ARM BOTTOM COVER & HARWARE PANTONE PQ COOL GRAY 1C ("MIST")		Yes	Yes	Yes	1	
KMT-031-000-201	MOUSE HOUSE		Yes	Yes	Yes	1	